

Product Datasheet

Goat Anti-Rat IgG H&L/Gold antibody (orb868899)

| | |
|------------------------|---|
| Description | Goat Anti-Rat IgG H&L/Gold antibody |
| Species/Host | Goat |
| Reactivity | Rat |
| Conjugation | Gold |
| Immunogen | Native Rat IgG |
| Preservatives | 10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol. |
| Form/Appearance | Liquid |
| Concentration | 0.4 mg/ml |
| Storage | Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Expiration Date | 12 months from date of receipt. |