

Product Datasheet

Rabbit Anti-Goat IgG H&L antibody (Gold) (orb868901)

Description	Rabbit Anti-Goat IgG H&L antibody (Gold)
Species/Host	Rabbit
Reactivity	Goat
Conjugation	Gold
Immunogen	Native Goat IgG
Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Form/Appearance	Liquid
Concentration	0.4 mg/ml
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.