

## Product Datasheet

### Rabbit Anti-Goat IgM antibody (Gold) (orb868910)

<b>Description</b>	Rabbit Anti-Goat IgM antibody (Gold)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Goat
<b>Conjugation</b>	Gold
<b>Immunogen</b>	Native Goat IgM
<b>Preservatives</b>	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.4 mg/ml
<b>Storage</b>	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.