

Product Datasheet

Mouse Anti-Rabbit IgM antibody (Gold) (orb868911)

| | |
|------------------------|-------------------------------------------------------------------------|
| Description | Mouse Anti-Rabbit IgM antibody (Gold) |
| Species/Host | Mouse |
| Reactivity | Rabbit |
| Conjugation | Gold |
| Immunogen | Native Rabbit IgM |
| Preservatives | 10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol. |
| Form/Appearance | Liquid |
| Concentration | 0.4 mg/ml |
| Storage | Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Expiration Date | 12 months from date of receipt. |