

Product Datasheet

C/EBPbeta (phospho-Thr235) antibody (FITC) (orb8690)

Description	C/EBPbeta (phospho-Thr235) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CEBP beta around the phosphorylation site of Thr235
Target	Phospho-CEBP beta (Thr235)
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P17676
Dilution Range	Flow-Cyt=1µg /test
Expiration Date	12 months from date of receipt.