

Product Datasheet

Rabbit Anti-Mouse IgM antibody (orb869151)

Description	Rabbit Anti-Mouse IgM antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	Native Mouse IgM
Preservatives	0.01M PBS (pH7.4)
Form/Appearance	Lyophilized powder
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB=1:1000-5000,IHC-P=1:100-500,IHC-F=1:200-1000,ELISA=1:2000-10000
Expiration Date	12 months from date of receipt.