

Product Datasheet

Mouse Anti-Rabbit IgM antibody (orb869152)

Description	Mouse Anti-Rabbit IgM antibody
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	Native Rabbit IgM
Preservatives	0.01M PBS (pH7.4)
Form/Appearance	Lyophilized powder
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC-P=1:100-500,IHC-F=1:200-1000,ELISA=1:2000-10000
Expiration Date	12 months from date of receipt.