

Product Datasheet

Mouse Anti-Rabbit IgG H&L ascites antibody (orb869159)

Description	Mouse Anti-Rabbit IgG H&L ascites antibody
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	Native Rabbit IgG
Form/Appearance	Lyophilized powder
Storage	Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB=1:10000-50000,IHC-P=1:100-500,IHC-F=1:10000-50000,ELISA=1:10000-100000
Expiration Date	12 months from date of receipt.