

Product Datasheet

Mouse Anti-Goat IgG H&L ascites antibody (orb869160)

Description	Mouse Anti-Goat IgG H&L ascites antibody
Species/Host	Mouse
Reactivity	Goat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	Native Goat IgG
Form/Appearance	Lyophilized powder
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB=1:10000-50000,IHC-P=1:100-500,IHC-F=1:10000-50000,ELISA=1:10000-100000
Expiration Date	12 months from date of receipt.