

## Product Datasheet

# Goat Anti-Rabbit IgG H&L/Streptavidin antibody (orb869551)

<b>Description</b>	Goat Anti-Rabbit IgG H&L/Streptavidin antibody
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	Streptavidin
<b>Tested Applications</b>	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
<b>Immunogen</b>	Native Rabbit IgG
<b>Preservatives</b>	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1.0 mg/ml
<b>Storage</b>	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	WB=1:500-2000,IHC-P=1:100- 500,IHC-F=1:100-500,ICC=1:100- 500,IF=1:100-500,ELISA=1:500-1000
<b>Expiration Date</b>	12 months from date of receipt.