

Product Datasheet

Mouse Anti-Goat IgG antibody (orb869555)

Description	Mouse Anti-Goat IgG antibody
Species/Host	Mouse
Reactivity	Goat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	Native Goat IgG
Preservatives	0.01M PBS (pH7.4)
Form/Appearance	Lyophilized powder
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,IF=1:100-500,ELISA=1:500-1000
Expiration Date	12 months from date of receipt.