

---

## Product Datasheet

**HER3 (phospho-Tyr1197) antibody (FITC)  
(orb8696)**

**Description**

Rabbit polyclonal to HER3 (phospho-

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	IF, IHC-Fr, IHC-P
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human HER3 around the phosphorylation site of Tyr1289 EE(p-Y)EY
<b>Target</b>	ERBB3
<b>Preservatives</b>	Lyophilized or Liquid
<b>Form/Appearance</b>	Lyophilized or Liquid
<b>Concentration</b>	- 100 µg (in 100 µl): 1 mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	148kDa
<b>CAS Number</b>	3
<b>Uniprot ID</b>	<b>P21860</b>
<b>Dilution Range</b>	IF=1:100-500, IHC-F=1:100-500, IHC-P=1:100-500
<b>Expiration Date</b>	12 months from date of receipt.