

Product Datasheet

RPA2 Antibody (PE/Cy5) (orb878786)

Description	Anti-phospho-RPA2(Thr21)/PE-Cy5
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	PE/Cy5
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RPA2 around the phosphorylation site of Thr21 GY(p-T)Q
Target	RPA2
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29 kDa
Uniprot ID	P15927
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.