

Product Datasheet

Phospho-PI3 Kinase p110 beta (Ser1070) Rabbit Polyclonal Antibody (Cy5) (orb881332)

Description	Phospho-PI3 Kinase p110 beta (Ser1070) Rabbit Polyclonal Antibody (Cy5)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy5
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PI3 Kinase p110 beta around the phosphorylation site of Ser1070 YR(p-S)
Target	PIK3CB
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ICC/IF=1:50-200, IF=1:100-500
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)