

Product Datasheet

APP Antibody (PE/Cy5) (orb889404)

Description	Anti-phospho-Amyloid Precursor Protein (Thr668)/PE-Cy5
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE/Cy5
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human APP around the phosphorylation site of Thr668 AV(p-T)PE
Target	APP
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	3.4-9/83kDa
Uniprot ID	P05067
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.