

Product Datasheet

phospho-NBN (Ser343)/AP Antibody (orb889445)

Description	Anti-phospho-p95 NBS1 (Ser343)/AP
Species/Host	Rabbit
Reactivity	Human
Conjugation	AP
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human NBS1 around the phosphorylation site of Ser343 SL(p-S)QG
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85kDa
Uniprot ID	O60934
Expiration Date	12 months from date of receipt.