

## Product Datasheet

# Glucocorticoid Receptor (phospho-Ser211) antibody (FITC) (orb8958)

<b>Description</b>	Glucocorticoid Receptor (phospho-Ser211) antibody (FITC)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC, IF
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human Glucocorticoid Receptor around the phosphorylation site of Ser211
<b>Target</b>	Phospho-Glucocorticoid Receptor (Ser211)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P04150</b>
<b>Dilution Range</b>	Flow-Cyt=1ug/Test,IF=1:100-500
<b>Expiration Date</b>	12 months from date of receipt.