
Product Datasheet

**Phospho-PTEN(Ser380 + Thr382 + Thr383)/PE-Cy3
antibody (orb900371)**

Description

Phospho-PTEN(Ser380 + Thr382 +

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy3
Tested Applications	FC, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PTEN around the phosphorylation site of Ser380/Thr382/Thr383
Target	Phospho-PTEN (Ser380 + Thr382 + Thr383)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44kDa
Uniprot ID	P60484
Dilution Range	Flow-Cyt=1µg /Test,ICC=1:100
Expiration Date	12 months from date of receipt.