

Product Datasheet

Human Cyclophilin D his tag Protein (orb90063)

Description	Human recombinant Cyclophilin-D, His Tag protein
Preservatives	1mg/ml solution containing 1x PBS pH-7.4 & 10% glycerol.
Form/Appearance	Sterile filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Enzymes
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MSHPSQAKP SNPSNPRVFFDVDIGGERVG RIVLELFADI VPKTAENFRA LCTGEKGIGH TTGKPLHFKG CPFHRIKKF MIQGGDFSNQ NGTGGESIYG EKFEDENFHY KHDREGLLSMANAGRNTNGS QFFITTVPTP HLDGKHVVFQ QVIKGIGVAR ILENVEVKGEKPAKLCVIAE CGELKEGDDG GIFPKDGSGD SHPDFPEDAD IDLKDVKIL LITEDLKNIG NTFKSNWE MAIKKYAEVL RYVDSSKAVI ETADRAKLQPIALSCVLNIG ACKLKMSNWQ GAIDSCLEAL ELDPSNTKAL YRRAQGWQGLKEYDQALADL KKAQGIAPED KAIQAELLKV KQKIKAQKDK EKAVYAKMFA
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.