

Product Datasheet

DENV Antibody (PE/Cy5) (orb900703)

Description	Anti-Dengue virus envelope glycoprotein E/PE-Cy5
Species/Host	Rabbit
Conjugation	PE/Cy5
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from Dengue virus type-2 envelope glycoprotein E (501-600/3390aa)
Target	DENV
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54kDa
Uniprot ID	P30026
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.