

Product Datasheet

Phospho-WNK1 (Thr60)/AP Antibody (orb901662)

Description	Anti-Phospho-WNK1(Thr60)/AP
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	AP
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human WNK1 around the phosphorylation site of Thr60 RH(p-T)MD
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	251kDa
Uniprot ID	Q9H4A3
Expiration Date	12 months from date of receipt.