

Product Datasheet

RPS6KA5 Antibody (PE/Cy7) (orb902982)

Description	Anti-Phospho-MSK1 (Ser360)/PE-Cy7
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	PE/Cy7
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MSK1 around the phosphorylation site of Ser360 TY(p-S)PA
Target	RPS6KA5
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90kDa
Uniprot ID	O75582
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.