

## Product Datasheet

### RPS6KA5 Antibody (PE/Cy5.5) (orb903309)

<b>Description</b>	Anti-Phospho-MSK1 (Ser376)/PE-Cy5.5
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	PE/Cy5.5
<b>Tested Applications</b>	FC
<b>Immunogen</b>	KLH conjugated synthesised phosphopeptide derived from human MSK1 around the phosphorylation site of Ser376 GY(p-S)FV
<b>Target</b>	RPS6KA5
<b>Preservatives</b>	Lyophilized or Liquid
<b>Form/Appearance</b>	Lyophilized or Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	90kDa
<b>Uniprot ID</b>	<a href="#">O75582</a>
<b>Dilution Range</b>	Flow-Cyt=2ug/Test
<b>Expiration Date</b>	12 months from date of receipt.