

Product Datasheet

RPS6KA5 Antibody (PE/Cy5.5) (orb903312)

Description	Anti-MSK1/PE-Cy5.5
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE/Cy5.5
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human MSK1 (31-130/802aa)
Target	RPS6KA5
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90kDa
Uniprot ID	O75582
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.