

Product Datasheet

Phospho-IRF3 (Ser396)/PE-Cy5.5 antibody (orb904176)

Description	Phospho-IRF3 (Ser396)/PE-Cy5.5 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	PE/Cy5.5
Tested Applications	ICC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human IRF3 around the phosphorylation site of Ser396
Target	Phospho-IRF3 (Ser396)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47kDa
Uniprot ID	Q14653
Expiration Date	12 months from date of receipt.