

## Product Datasheet

### APP Antibody (PE/Cy5.5) (orb905413)

<b>Description</b>	Anti-phospho-Amyloid Precursor Protein (Thr668)/PE-Cy5.5
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	PE/Cy5.5
<b>Tested Applications</b>	WB
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human APP around the phosphorylation site of Thr668 AV(p-T)PE
<b>Target</b>	APP
<b>Preservatives</b>	Lyophilized or Liquid
<b>Form/Appearance</b>	Lyophilized or Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	3.4-9/83kDa
<b>Uniprot ID</b>	<b>P05067</b>
<b>Dilution Range</b>	WB=1:500-2000
<b>Expiration Date</b>	12 months from date of receipt.