

## Product Datasheet

# Laemmli Buffer: SDS-PAGE Sample Loading Buffer 2X (Reducing) (orb90557)

<b>Description</b>	Laemmli Buffer: SDS-PAGE Sample Loading Buffer 2X (Reducing)
<b>Storage</b>	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.. Shipped in room temperature.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.