

## Product Datasheet

### APP Antibody (PE/Cy5) (orb906543)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Anti-beta-Amyloid (1-42)/PE-Cy5   |
| <b>Species/Host</b>        | Mouse   |
| <b>Conjugation</b>         | PE/Cy5  |
| <b>Tested Applications</b> | WB  |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human beta-Amyloid (1-42/42aa)    |
| <b>Target</b>              | APP   |
| <b>Preservatives</b>       | Lyophilized or Liquid   |
| <b>Form/Appearance</b>     | Lyophilized or Liquid   |
| <b>Concentration</b>       | 1mg/ml  |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | 4.4kDa  |
| <b>Uniprot ID</b>          | <a href="#">P05067</a>  |
| <b>Dilution Range</b>      | WB=1:500-2000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |