

Product Datasheet

phospho-AKT3 (Tyr473)/PE-Cy5 Antibody (orb908305)

Description	Anti-phospho-AKT3 (Tyr473)/PE-Cy5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy5
Tested Applications	FC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human AKT3 around the phosphorylation site of Tyr473 FS(p-Y)SA
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57kDa
Uniprot ID	P31749
Expiration Date	12 months from date of receipt.