

Product Datasheet

APP Antibody (PE/Cy5) (orb909387)

Description	Anti-beta-Amyloid (1-42)/PE-Cy5
Species/Host	Mouse
Reactivity	Rat
Conjugation	PE/Cy5
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human beta-Amyloid (672-713/770 or 1-42/42aa)
Target	APP
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	4.4kDa
Uniprot ID	P05067
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.