

Product Datasheet

APP Antibody (PE/Cy5.5) (orb912517)

Description	Anti-beta-Amyloid 1-40 (CT)/PE-Cy5.5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy5.5
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human beta-Amyloid 1-40 (37-40/40aa)
Target	APP
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	4.3kDa
Uniprot ID	P05067
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.