

Product Datasheet

APP Antibody (PE/Cy5.5) (orb913049)

Description	Anti-beta Amyloid 1-16/PE-Cy5.5
Species/Host	Mouse
Reactivity	Human, Rabbit, Rat
Conjugation	PE/Cy5.5
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide of human beta-Amyloid (1-16/42aa)
Target	APP
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	4.3kDa
Uniprot ID	P05067
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.