
Product Datasheet

LKB1 (phospho-Ser428) antibody (FITC) (orb9170)

Description

Rabbit polyclonal to LKB1 (phospho-

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human LKB1 around the phosphorylation site of Ser428 [RL(p-S)AC]
Target	STK11
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48kDa
Uniprot ID	Q15831
Dilution Range	IF: 1:50-200
Expiration Date	12 months from date of receipt.