

Product Datasheet

APP Antibody (PE/Cy7) (orb917310)

Description	Anti-Amyloid Precursor Protein/PE-Cy7
Species/Host	Mouse
Reactivity	Rat
Conjugation	PE/Cy7
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from the middle of human Soluble APP-beta/alpha (381-480/770aa)
Target	APP
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	72-83kDa
Uniprot ID	P05067
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.