

Product Datasheet

MAML3/AP Antibody (orb920067)

Description	Anti-MAML3/AP
Species/Host	Rabbit
Reactivity	Human
Conjugation	AP
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human MAML3 (931-1030/1134aa)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	122kDa
Uniprot ID	Q96JK9
Expiration Date	12 months from date of receipt.