

Product Datasheet

phospho-C-REL (Ser503)/PE-Cy7 Antibody (orb923341)

Description	Anti-phospho-C-REL (Ser503)/PE-Cy7
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	ICC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human C-REL around the phosphorylation site of Ser503 TS(p-S)DS
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	69kDa
Uniprot ID	Q04864
Expiration Date	12 months from date of receipt.