

Product Datasheet

Nrf2 (phospho-Ser40) antibody (FITC) (orb9392)

Description	Nrf2 (phospho-Ser40) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	FITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Nrf2 around the phosphorylation site of Ser40
Target	phospho-Nrf2 (Ser40)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16236
Dilution Range	Flow-Cyt=1µg /test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.