

## Product Datasheet

### Random Hexamers (orb93926)

<b>Description</b>	Random Hexamers are a mixture of oligonucleotides representing all possible sequences for a hexamer...
<b>Form/Appearance</b>	clear aqueous solution
<b>Concentration</b>	100 $\mu$ M
<b>Storage</b>	store at -20 °C
<b>Note</b>	For research use only
<b>Purity</b>	by HPLC
<b>MW</b>	Theoretical MW: 1793 g/mol
<b>Expiration Date</b>	12 months from date of receipt.