

Product Datasheet

DENV Antibody (PE/Cy7) (orb942590)

Description	Anti-Dengue Virus NS1 /PE-Cy7
Species/Host	Rabbit
Conjugation	PE/Cy7
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human DENV-2 NS1 (801-900/3391aa)
Target	DENV
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40kDa
Uniprot ID	P30026
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.