

---

## Product Datasheet

**p70 S6 Kinase Beta (phospho-Thr389) antibody  
(FITC) (orb9465)**

**Description**

Rabbit polyclonal to p70 S6 Kinase

conjugated to FITC

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse

**Conjugation**

FITC

**Tested**

FC

**Applications****Immunogen**

KLH conjugated Synthesised phosphopeptide derived from human p70 S6 Kinase Beta 2 around the phosphorylation site of Tyr389

**Target**

RPS6KB2

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

- 100 µg (in 100 µl): 1 mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

53kDa

**Uniprot ID**[Q9UBS0](#)**Dilution Range**

IF: 1:50-200

**Expiration Date**

12 months from date of receipt.