


Product Datasheet

Human CXCL5/LIX/ENA-78 ELISA Kit (orb96782)

Description	Human CXCL5/LIX/ENA-78 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A37-N114
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	For serum, dilution ratio is 1:1 and concentration is 442-3544 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 0-1082 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	11972 MW
Uniprot ID	P42830
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human CXCL5/ENA-78 PicoKine ELISA Kit st...