


## Product Datasheet

### Human HBEGF ELISA Kit (orb96788)

<b>Description</b>	Human HBEGF ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	7.8 pg/ml - 500 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: NS0; Immunogen sequence: D63-L148
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Application notes</b>	For serum, dilution ratio is 1:1 and concentration is 36-451 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 0-114 pg/ml
<b>Sample Types</b>	culture supernatants, cell lysates, serum, plasma (heparin, EDTA), saliva, urine and human milk.
<b>MW</b>	23067 MW
<b>Uniprot ID</b>	<a href="#">Q99075</a>
<b>Sensitivity</b>	10 pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
HBEGF  
PicoKine  
ELISA Kit  
standard ...