


Product Datasheet

Human Clusterin ELISA Kit (orb96800)

Description	Human Clusterin ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D23-R227 (beta) & S228-E449 (alpha)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:50000, concentration in serum and plasma is 156-453 ug/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva and urine.
MW	52495 MW
Uniprot ID	P10909
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Human Clusterin PicoKine ELISA Kit stand...