


## Product Datasheet

# Mouse CXCL5 / LIX / ENA-78 / GCP 2 ELISA Kit (orb96802)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Mouse CXCL5 / LIX / ENA-78 / GCP 2 ELISA Kit   |
| <b>Reactivity</b>          | Mouse  |
| <b>Range</b>               | 15.6 pg/ml - 1,000 pg/ml   |
| <b>Tested Applications</b> | ELISA  |
| <b>Immunogen</b>           | Expression system for standard: E.coli;<br>Immunogen sequence: V45-A118  |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months.<br>Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only  |
| <b>Sample Types</b>        | cell culture supernatants, serum and plasma (heparin, EDTA, citrate).  |
| <b>MW</b>                  | 14235 MW   |
| <b>Uniprot ID</b>          | <b>P50228</b>  |
| <b>Sensitivity</b>         | 10 pg/ml   |
| <b>Expiration Date</b>     | Please enquire.  |



Mouse CXCL5/ENA-78 PicoKine ELISA Kit st...