


Product Datasheet

Mouse Clusterin / Apolipoprotein J ELISA Kit (orb96806)

Description	Mouse Clusterin / Apolipoprotein J ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D22-R226 (beta) & S227-E448 (alpha)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:50000, concentration in serum and plasma is around 100 ug/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	51656 MW
Uniprot ID	Q06890
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse Clusterin PicoKine ELISA Kit stand...