


Product Datasheet

Mouse DKK-1 ELISA Kit (orb96807)

Description	Mouse DKK-1 ELISA Kit
Reactivity	Mouse
Range	46.9 pg/ml - 3,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: S30-H272
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:20, concentration in serum and plasma is 1156-33076 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	29298 MW
Uniprot ID	054908
Sensitivity	15 pg/ml
Expiration Date	Please enquire.



Mouse DKK1
PicoKine
ELISA Kit
standard c...