

---

## Product Datasheet

### BCRP/ABCG2 Antibody (orb96996)

**Description**

BCRP/ABCG2 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

WB

**Immunogen**

A synthetic peptide corresponding to a sequence in the middle region of human BCRP/ABCG2 (154-168aa NHEKNERINRVIQEL).

**Form/Appearance**

Lyophilized

**Concentration**

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

**Storage**

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Application notes**

WB: The detection limit for ABCG2 is approximately 1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

**Isotype**

Rabbit IgG

**Clonality**


Polyclonal

**MW**

72314 MW

**Uniprot ID**
[Q9UNQ0](#)
**Expiration Date**

12 months from date of receipt.



Western blot analysis of Lane 1: HELA Ce...