

## Product Datasheet

### BRMS1/PE-Cy5 Antibody (orb973077)

<b>Description</b>	Anti-BRMS1/PE-Cy5
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	PE/Cy5
<b>Tested Applications</b>	ICC
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from human BRMS-1 (151-246/246aa)
<b>Preservatives</b>	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Lyophilized or Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	28kDa
<b>Uniprot ID</b>	<a href="#">Q9HCU9</a>
<b>Expiration Date</b>	12 months from date of receipt.