



Product Datasheet

SERPINE1 Antibody (orb97411)


| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to SERPINE1 |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FC, IHC |
| Immunogen | Purified recombinant fragment of human SERPINE1 expressed in E. Coli. |
| Preservatives | Purified antibody in PBS with 0.05% sodium azide |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | Mouse IgG1 |
| Clonality | Monoclonal |
| MW | 45kDa |
| Entrez | 5054 |
| Dilution Range | IHC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000 |
| Expiration Date | 12 months from date of receipt. |



Flow cytometric analysis of NIH/3T3 cell...



IHC-P of kidney cancer tissues using SER...



Western blot analysis of HEK293 (Lane1) ...