



Product Datasheet

PPP1CA Antibody (orb97412)


Description	Mouse monoclonal antibody to PPP1CA
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant fragment of human PPP1CA expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
IsoType	Mouse IgG1
Clonality	Monoclonal
MW	38kDa
Entrez	5499
Dilution Range	WB: 1/500 - 1/2000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human PPP1CA (A...



Line graph illustrates about the Ag-Ab r...



Western blot analysis of Hela (Lane1), H...